

DESCOCEPT PUR



Alcoholic rapid
disinfectant



PRESERVE

PRODUCT DESCRIPTION

- Ready to use
- Residue-free
- Fragrance-free

DESCOCEPT PUR is a ready-to-use, aldehyde-free rapid disinfectant based on ethanol for a rapid and effective disinfection of medical inventory, medical devices as well as all kinds of surfaces. Especially suitable for difficult to reach surfaces.

DESCOCEPT PUR is free from perfume and residue-free. It is therefore especially applicable for the rapid disinfection in the food industry and in canteen kitchens..

RANGE OF APPLICATION

According to Biocidal Products Regulation (BPR)

Disinfection and cleaning of alcohol resistant medical inventory and all kinds of surfaces. Due to the perfume-free and residue-free formulation ideally suited for use in food industry and canteen kitchens.

Acc. to EU Medical Devices Regulation

Disinfection and cleaning of non-invasive alcohol-resistant medical devices.

Application

To ensure complete wetting apply the undiluted solution evenly on alcohol-resistant surfaces and keep them moist for the required contact time.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

DESCOCEPT PUR can be used for wipe disinfection or as a spray.

Use in dispenser systems: DESCOCEPT PUR is especially suited for the use with wipe dispenser systems for a convenient and user-friendly application with a proven disinfectant efficacy and standing time of up to 42 days in combination with the ONE SYSTEM PLUS wipe dispensers and with a standing time up to 28 days in combination with the DESCO/ECO WIPES dispenser system.

Composition

100 g solution contain: 45 g ethanol.

Material compatibility

The use is not recommended for the following materials: polymethylmethacrylate (PMMA) and polyethylene terephthalate (PET).

For application on the following materials we recommend only under certain conditions: copper, brass and acrylonitrile-butadiene-styrene (ABS).

Precautionary and hazard statements

Flammable liquid and vapour. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Do not breathe vapour/spray. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents/container to approved disposal company or local collection.

For professional use only by personnel with corresponding specialist knowledge according to national directives.

Read label before use.

DESCOCEPT PUR



PRESERVE

SPECTRUM OF ACTIVITY			30 s	1 min	3 min	5 min
Application recommendation for surface disinfection						
bactericidal, yeasticidal*	VAH / EN 16615	with mechanical action, dirty condition				•
	VAH / EN 13697	without mechanical action, clean and dirty conditions				•
tuberculocidal (M. terrae)	EN 14348	clean condition	•			
mycobactericidal (M. avium and M. terrae)	EN 14348	clean condition		•		
limited spectrum virucidal activity	EN 14476	clean condition			•	
	EN 16777	clean condition, without mechanical action				•
active against enveloped viruses	EN 14476	clean and dirty conditions		•		
	EN 16777	clean and dirty conditions, without mechanical action		•		
Additional test results						
bactericidal	EN 13727	clean and dirty conditions		•		
	EN 1276	clean and dirty conditions		•		
	EN 13697	clean and dirty conditions, without mechanical action		•		
yeasticidal (C. albicans)	EN 13624	clean and dirty conditions		•		
	EN 1650	clean and dirty conditions		•		
	EN 13697	clean and dirty conditions, without mechanical action		•		
activity against enveloped viruses	RKI / DVV	incl. HBV/HIV/HCV	•			
active against rotaviruses	EN 14476	clean condition	•			
active against adenoviruses	EN 14476	clean condition			•	
	EN 16777	clean condition, without mechanical action				•
active against noroviruses (MNV)	EN 14476	clean condition		•		
	EN 16777	clean condition, without mechanical action				•

*einschließlich Phase 2 Stufe 1 und Phase 2 Stufe 2 Tests (quantitative Suspensionsversuche und praxisnahe Keimträgerversuche)

Product	Single unit	Delivery unit	Content	Art. Nr.
DESCOCEPT PUR	Bottle	10	1L	00-323-010E
	Canister	3	5L	00-323-050E

National specifications may differ. For further information, please contact our subsidiary or your local dealer.
Use disinfectants safely. Always read label and product information before use.



CE 0123

CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:

<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001, has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.