

CLEANISEPT®



PREVENT

SURFACE DISINFECTION AND CLEANING FOR NON-INVASIVE MEDICAL DEVICES WITH MECHANICAL ACTION

1 / 9

Good disinfection & cleaning performance

Convinces with a good material compatibility
even on alcohol-sensitive surfaces

No addition of
colorants and fragrances



SURFACE / DISINFECTANT CONCENTRATES

CLEANISEPT®



PREVENT

PRODUCT DESCRIPTION

2 / 9

CLEANISEPT is a surface disinfectant and cleaner for non-invasive medical devices with mechanical action. **CLEANISEPT** is based on quaternary ammonium compounds. It is fragrance-free, contains neither aldehydes nor phenols and is characterized by good cleaning strength.

CLEANISEPT is suitable for surface disinfection of washable non-porous surfaces and for areas in which odor nuisance must be avoided, in addition to meeting high hygienic standards. **CLEANISEPT** has good material compatibility and is also suitable for acrylic glass.

APPLICATIONS AND NOTES

Range of application

Surface disinfection and cleaning for non-invasive medical devices with mechanical action.

According to Biocidal Products Regulation (BPR): For surface disinfection and cleaning of medical inventory and all kinds of surfaces. Suitable for use in the food sector.

Application

For professional use only. Visible contamination must be removed before disinfection. Wet the surface to be disinfected completely with a sufficient amount of the freshly prepared **CLEANISEPT** working solution and spread it mechanically (scrub-wipe disinfection).

Concentrate - only use diluted! Prepare the desired concentration of the usage solution with water of suitable quality (max. 40 °C).

Medical devices that come into contact with the body during intended use should be wiped down with a wet wipe (water of suitable quality) after the exposure time has elapsed.

Thoroughly rinse surfaces and tools that come into contact with food with clean water (drinking water quality) after the exposure time.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Use with wipe dispenser systems: **CLEANISEPT** is well suited for the use in wipe dispenser systems with a proven disinfectant efficacy and standing time of up to 42 days in combination with the ONE SYSTEM wipe dispensers (PET non-woven) or with a standing time of up to 28 days in combination with the DESCO/ECO Wipes dispenser system (PET non-woven).

Composition

100 g contain: 6.66 g alkyl (C12-16) dimethylbenzylammonium chloride, 3.33 g didecylmethylammonium chloride.

Material compatibility

Wide range of applications on surfaces and medical devices. (see page 4 – 7)

Precautionary and hazard statements

Causes skin irritation. Causes serious eye damage. Very toxic to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Dispose of contents/container to an approved waste disposal plant.

Use disinfectants safely.

Always read label and product information before use.

SURFACE / DISINFECTANT CONCENTRATES

CLEANISEPT®



PREVENT

SPECTRUM OF EFFICACY AND CONTACT TIMES

3 / 9



RECOMMENDATION FOR SURFACE DISINFECTION WITH MECHANICAL ACTION				1 min	5 min	15 min	30 min	60 min
application recommendation for surface disinfection								
bactericidal, yeasticidal	EN 13727 EN 13624 EN 16615	dirty conditions		7,5 %		2,5 %	1 %	
virucidal against enveloped viruses	EN 14476	clean and dirty conditions		7,5 %		2 %	1,25 %	
active against noroviruses	EN 14476	dirty conditions						7,5 %
Recommendation for surface disinfection with mechanical action in non-medical areas								
bactericidal	EN 1276	clean and dirty conditions			0,5 %			
yeasticidal	EN 1650	dirty conditions			0,5 %			

CLEANISEPT®



PREVENT

MATERIAL COMPATIBILITY

4 / 9



MATERIAL METALS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
stainless steel V2A			•	Medical transport chairs
			•	Toilet chairs
stainless steel V4A			•	Rolling walking aids
aluminum			•	

CLEANISEPT®



PREVENT

MATERIAL COMPATIBILITY

5 / 9



MATERIAL PLASTICS: ELASTOMERS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
PUR (polyurethane)			•	Medical transport chairs
silicones			•	Respiratory masks
			•	Medical keyboards and mice
Latex rubber			•*	
EPDM (ethylene propylene diene (monomer) rubber)			•	
CR (neoprene)			•	

* Test the material compatibility of these materials in advance on a small, inconspicuous area.

CLEANISEPT®



PREVENT

MATERIAL COMPATIBILITY

6 / 9



MATERIAL PLASTICS: THERMOPLASTICS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
PE-HD (polyethylene-high density)			•	Storage and transport containers
POM (polyoxymethylene)			•	Plastic chairs
PA (polyamide)			•	X-ray apron
PSU (polysulfone)			•	Housing components
PVC (polyvinylchloride)			•	
PMMA (polymethylmethacrylate)			•	Acrylic and plexiglass, incubators
PC (polycarbonate e.g. Makrolon)	•			Ultrasound devices
	•			ECG devices
ABS (acrylonitril-butadiene-styrene)	•			Patient monitoring monitors
	•			Medical keyboards and mice

SURFACE / DISINFECTANT CONCENTRATES

CLEANISEPT®



PREVENT

MATERIAL COMPATIBILITY

7 / 9



MATERIAL TEXTILES AND SURFACES	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
synthetic leather (e.g. with PVC and PUR top layer)			* Examination and treatment couches * Gynaecological examination chair	
			* Dental chairs	
			* Operating couches	
PVC (polyvinylchloride) (2 %)			Floor	
wood floor (2 %)				
linoleum (2 %)			* Floor	

* Test the material compatibility of these materials in advance on a small, inconspicuous area.

SURFACE / DISINFECTANT CONCENTRATES

CLEANISEPT®



PREVENT

MANUFACTURER APPROVALS

8 / 9



Honeywell

PURCHASING INFORMATION

pH 0 1 2 3 4 5 6 7 8 9 10 11 12 13 14  CE 0123

Product	Single Unit	Delivery Unit	Content	REF
CLEANISEPT	dosage bottle	10	1 L	00-204-010
	bottle	6	2 L	00-204-020
	canister	3	5 L	00-204-050

The availability of the products and container sizes depend on a successful national registration.
National specifications may differ. For further information, please contact our subsidiary or your local dealer.

CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:
<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.

SURFACE / DISINFECTANT CONCENTRATES

CLEANISEPT®



PREVENT

PRODUCT FAMILY OVERVIEW

9 / 9



CLEANISEPT WIPES



CLEANISEPT WIPES MAXI



CLEANISEPT WIPES



CLEANISEPT



RECOMMENDED NON-WOVEN WIPE DISPENSER SYSTEMS



ONE SYSTEM+ PLUS



ONE SYSTEM BASIC



DESCO WIPES



ECO WIPES