

SURFACE

DESCOCEPT PUR WIPES RTU



PRESERVE

ALCOHOL-IMPREGNATED WIPES FOR RAPID DISINFECTION

1 / 6

Versatile application, **comfortable and time saving**

Extensive, rapid efficacy

Especially suited **for disinfection of large-sized surfaces** in near patient areas

Residue-free rapid disinfection



SURFACE

DESCOCEPT PUR WIPES RTU



PRESERVE

PRODUCT DESCRIPTION

2 / 6

DESCOCEPT PUR WIPES RTU are alcoholic, ready to use wipes for disinfection and cleaning of non invasive medical devices and medical inventory as well as surfaces in patient-near areas.

DESCOCEPT PUR WIPES RTU offer a universal and time saving application for the disinfection and cleaning of large surfaces. The wipes can be extracted individually

from the known wipes dispenser system with re-sealable lid. Each refill pack comes with an ULTRASOL OXY SACHET to ensure proper, manual reprocessing before refilling the dispenser.

APPLICATIONS AND NOTES

Acc. to Biocidal Products Regulation (BPR)

Rapid disinfection and cleaning of alcohol-resistant medical inventory and all types of surfaces, mainly in near patient areas.

Acc. to EU Medical Devices Regulation

Rapid disinfection and cleaning of alcohol resistant, non-invasive medical devices.

Application

Wipe surfaces with DESCOCEPT PUR WIPES RTU until completely wet and keep them moist for the required contact time.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Shelf life after opening: 28 days.

Composition

100 g solution contain: 45 g ethanol.

Material compatibility

Wide range of applications on surfaces and medical devices. (see page 4 – 5)

Product status

Dual Registration (medical device/biocide)

Precautionary and hazard statements

Flammable liquid and vapour. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents/container to approved disposal company or local collection.

For professional use only by personnel with corresponding specialist knowledge according to national directives.

Use disinfectants safely.

Always read label and product information before use.

SURFACE

DESCOCEPT PUR WIPES RTU



PRESERVE

SPECTRUM OF EFFICACY AND CONTACT TIMES

3 / 6



SPECTRUM OF EFFICACY AND CONTACT TIMES			30 s	1 min	3 min	5 min
application recommendation for surface disinfection						
bactericidal, yeasticidal ¹	VAH/EN 16615	with mechanical action, dirty conditions				•
tuberculocidal (M. terrae)	EN 14348	clean conditions	•			
mycobactericidal (M. avium und M. terrae)	EN 14348	clean conditions		•		
limited spectrum virucidal activity	EN 14476	clean conditions			•	
active against enveloped viruses	EN 14476	clean and dirty conditions		•		
additional test results						
active against adenoviruses	EN 14476	clean conditions			•	
active against noroviruses (MNV)	EN 14476	clean conditions		•		
bactericidal	EN 13727	clean and dirty conditions		•		
	EN 1276	clean and dirty conditions		•		
yeasticidal (C. albicans)	EN 13624	clean and dirty conditions	•			
	EN 1650	clean and dirty conditions		•		

¹ – including phase 2 step 1 and phase 2 step 2 tests (quantitative suspension tests and quantitative carrier tests)

SURFACE

DESCOCEPT PUR WIPES RTU



PRESERVE

MATERIAL COMPATIBILITY

4 / 6



MATERIAL METALS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
Stainless steel V2A			•	Medical transport chairs
			•	Toilet chairs
Stainless steel V4A			•	Rolling walking aids
			•	Walking frames
Aluminum			•	
Copper		•		
Brass		•		Luer lock adapter
		•		Stethoscope (chrome-plated brass)

SURFACE

DESCOCEPT PUR WIPES RTU



PRESERVE

MATERIAL COMPATIBILITY

5 / 6



MATERIAL PLASTICS: THERMOPLASTICS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
PC (polycarbonate e.g. Makrolon)			• Ultrasound devices • ECG devices • EEG devices • MRI devices • CT devices	
POM (polyoxymethylene)			• Plastic chairs	
ABS (acrylonitril-butadiene-styrene)		• Patient monitoring monitors • Medical keyboards and mice		
PET (polyethylene terephthalate)	•			
PMMA (polymethylmethacrylate)	•		• Acrylic and plexiglass, incubators	

SURFACE

DESCOCEPT PUR WIPES RTU



PRESERVE

PURCHASING INFORMATION

6 / 6

Product	Single Unit	Delivery Unit	Content	Wipe Size	Art. no.
DESCOCEPT PUR WIPES RTU	refill	2	120 wipes	17,5 x 36 cm	00-323-RS120E

The availability of the products and container sizes depend on a successful national registration.

pH 0 1 2 3 4 5 6 7 8 9 10 11 12 13 14   0123

CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:

<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.