

INSTRUMENTS / MAUNAL REPROCESSING

DESCO CLEANER



INSTRUMENT CLEANER

1 / 3

Thorough and quick cleaning
Suitable for cleaning procedures using ultrasound
Contains biodegradable surfactants



INSTRUMENTS / MAUNAL REPROCESSING

DESCO CLEANER



PRODUCT DESCRIPTION

2 / 3

DESCO CLEANER is an alkaline powder cleaner for instruments for thorough quick cleaning in medical areas and laboratories. For heavily soiled surfaces with protein and blood residues, **DESCO CLEANER** ensures efficient and easy cleaning of surgical instruments and medical devices.

APPLICATIONS AND NOTES

Range of application

Intended purpose acc. to EU Medical Devices Regulation: Cleaner for instruments and medical devices.

Application

For professional use only. Suitable for alkaline resistant dental, medical and surgical instruments, as well as laboratory equipment.

Working solution in the desired concentration with water (max. 40 °C). Prepare solution by diluting with suitable water in accordance with national guidelines. Immerse open instruments and devices in the working solution directly after use and leave to soak for 10-20 min. Every surface and cavity to be cleaned must be wetted entirely by the working solution. In case of persistent soiling, increase the dose and/or the exposure time. If necessary, clean mechanically, e.g. by brushing. After cleaning, the instruments must be removed from the solution, thoroughly rinsed with suitable water and further processed. DESCO CLEANER is suitable for the use in ultrasound baths.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Prepare solution by diluting with suitable water (max. 40 °C) in accordance with national guidelines. For cleaning in the ultrasonic bath: Do not exceed the ultrasonic exposure time according to the manufacturer's instructions or temperatures of 40°C. Change working solution immediately if there is any visible contamination but at least once every working day. DESCO CLEANER must not be mixed with other process chemicals.

Dosage

Normal soiling: 1/2 measuring spoon into 3 litres of water (0.5 %)

High soiling: 1 measuring spoon into 3 litres of water (1 %)

Exposure time: 10-20 min; for ultrasound bath: 5-10 min.

Composition

Phosphates > 30 %, Non-ionic surfactants 5-15 %.

The surfactants contained in the product are biodegradable according to the EU Regulation on Detergents.

Material compatibility

The use of DESCO CLEANER is not recommended for the following materials: polycarbonate (PC).

Polymethylmethacrylate (PMMA), polyethylene terephthalate (PET), polysulfone (PSU), chloroprene (CR/SBR) and ethylene propylene diene (monomer) rubber (EPDM) may undergo superficial changes that do not affect functionality.

Precautionary and hazard statements

May be corrosive to metals. Causes severe skin burns and eye damage. May cause respiratory irritation. Do not breathe dust. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Store in a corrosion resistant container with a resistant inner liner. Dispose of contents/container to an approved waste disposal plant.

Read label before use.

INSTRUMENTS / MAUNAL REPROCESSING

DESCO CLEANER



PURCHASING INFORMATION

3 / 3

Product	Single unit	Delivery unit	Content	REF
DESCO CLEANER	Box	6	1 kg	00-115-010
	Bucket	1	5 kg	00-115-050

National information may differ. For further information, please contact our subsidiary or your local dealer.
The availability of the products and container sizes depend on a completed national registration.



CERTIFICATIONS



Unter folgendem Link finden Sie unsere Zertifikate und Siegel:

<https://www.schumacher-online.com/de/zert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.