

MANUAL CLEANING

ORANGE SOLVENT



CLEANER FOR SURFACE CLEANING OF MEDICAL DEVICES

1 / 3

With natural orange oil
With extra cleaning power
Effortlessly removes sticky residues and markings



MANUAL CLEANING

ORANGE SOLVENT



PRODUCT DESCRIPTION

2 / 3

ORANGE SOLVENT is a ready-to-use cleaning solution for thorough removal of stubborn residues on medical instruments and equipment.

Also for removal of adhesive residue, permanent marker and scuff marks on flooring. Economic to use.

APPLICATIONS AND NOTES

Range of application

Cleaner for surface cleaning of medical devices.

Application

For professional use only. Apply ORANGE SOLVENT to a wipe or single-use cloth and rub the area to be treated.

When used on medical devices, these must subsequently be fully reprocessed in a validated process in accordance with national regulations.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

ORANGE SOLVENT must not be mixed with other process chemicals.

Easy removal of stubborn stickers and old labels after a short application time and easy removal of adhesive residue.

- For substantial removal of iodine specks and other colour signal specks.
- For removing stubborn scuff marks on PVC and linoleum flooring such as those caused by bed castors.

Composition

Citrus Sinensis Peel Extract, Limonene, Decyl Oleate.

Material compatibility

The use is not recommended for the following materials: polycarbonate (PC), polymethylmethacrylate (PMMA), acrylonitrile butadiene styrene (ABS), polyoxymethylene (POM)

Do not use with disposable plastic cups and disposable latex gloves.

Precautionary and hazard statements

Flammable liquid and vapour. May be fatal if swallowed and enters airways. Causes skin irritation. May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Immediately call a POISON CENTER/doctor. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. Do NOT induce vomiting. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents/container to an approved waste disposal plant.

Read label before use.

MANUAL CLEANING

ORANGE SOLVENT



PURCHASING INFORMATION

3 / 3

Product	Single unit	Delivery unit	Content	REF
ORANGE SOLVENT	bottle	9	250 mL	00-805-0025
	bottle	6	500 mL	00-805-005

National information may differ. For further information, please contact our subsidiary or your local dealer.
The availability of the products and container sizes depend on a completed national registration.



CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:
<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.