

PERFEKTAN® TB



PREVENT

ALDEHYDE-FREE CONCENTRATE FOR INSTRUMENT CLEANING AND DISINFECTION

1 / 4

Aldehyde-free preparation with good cleaning power
Gentle on materials
Pleasant scent



PERFEKTAN® TB



PRODUCT DESCRIPTION

2 / 4

PERFEKTAN TB is an economical, material-friendly, cleaning instrument disinfectant for routine clinical use. Thanks to the antimicrobial efficacy personal protection is assured supported instrument processing. **PERFEKTAN TB** is characterised by a pleasant smell thanks to the aldehyde-free formulation and is suitable for use on medical instruments.

APPLICATIONS AND NOTES

Range of application

For non-fixing disinfecting cleaning and for disinfection of instruments and dental instruments. Suitable for use with ultrasonic baths. Not suited for reprocessing of laryngeal masks.

Application

Prepare **PERFEKTAN TB** working solution in the desired concentration with water (max. 40 °C). Prepare solution by diluting with suitable water in accordance with national guidelines.

Disinfectant cleaning in the immersion bath procedure: Place instruments, after any preliminary cleaning that may be necessary, into the working solution and clean mechanically where required. After the exposure time has elapsed, rinse the instruments thoroughly with water and continue processing. Suitable for use in ultrasonic baths.

Disinfection in the immersion bath procedure: Clean the instruments and place them into the working solution. All surfaces and cavities must be wetted completely by the working solution. After the exposure time has elapsed, rinse the instruments thoroughly with suitable water, dry and continue processing.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

For cleaning in the ultrasonic bath: Do not exceed the sonication time as specified in the manufacturer's instructions or temperatures of 40 °C. Not suitable for final disinfection of semi-critical and critical medical devices. Change working solution if there is any visible contamination but at least once every working day. **PERFEKTAN TB** must not be mixed with other process chemicals. For final disinfection of medical devices with a risk class of semi-critical A or higher, process using a virucidal disinfectant, e.g. **PERFEKTAN ACTIVE**. Ultrasonic bath procedure: **PERFEKTAN TB** can be used to improve cleaning performance in all ultrasonic baths made from stainless steel. The freshly prepared **PERFEKTAN TB** working solution used for cleaning must be de-gassed before first use. Avoid temperatures greater than 40 °C in the working solution, because these could result in encrustation. After the contact

time has elapsed, rinse the medical devices thoroughly with demineralised water. Appropriate personal protective equipment (PPE) must be worn when handling **PERFEKTAN TB**. Before using **PERFEKTAN TB** for the first time, clean all instruments and baths thoroughly to remove chemical residues. Instrument baths must be cleaned and disinfected at least once every working day, preferably with a disposable wipe soaked in alcohol, e.g. **DESCOCEPT PUR WIPES RTU**.

Composition

100 g contains: 5.6 g N, N-didecyl-N-methyl-poly(oxyethyl)ammonium propionate; 3,75 g cocospropylenediamine-1,5-bis-guanidinium diacetate. Contains fragrances: Limonene, 5 – 15 % non-ionic surfactants, disinfectant.

Material compatibility

The use is not recommended for the following materials: brass, aluminum, polystyrene (PS), polysulfone (PSU), polyethylen terephthalate (PET), polymethylmethacrylate (PMMA), acrylonitrile-butadiene-styrene (ABS), polycarbonate (PC). For application on the following materials we recommend only under certain conditions: copper, stainless steel V2A.

Precautionary and hazard statements

Causes skin irritation. Causes serious eye damage. Very toxic to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: Call a POISON CENTER/doctor. Dispose of contents/container to approved disposal company or local collection.

For professional use only by personnel with corresponding specialist knowledge according to national directives.

Read label before use.

PERFEKTAN® TB



SPECTRUM OF EFFICACY AND CONTACT TIMES

3 / 4

SPECTRUM OF EFFICACY AND CONTACT TIMES			5 min	15 min	30 min	60 min
APPLICATION RECOMMENDATION						
bactericidal, yeasticidal	EN 14561 EN 14562	dirty conditions		3 %	2 %	1 %
active against enveloped viruses	EN 14476 EN 17111	dirty conditions		3 %		
ADDITIONAL TEST RESULTS						
bactericidal ¹ , yeasticidal ¹	VAH/EN ²	dirty conditions			2 %	1 %
bactericidal	EN 13727	dirty conditions	1,5 %	1 %		
	EN 14561	dirty conditions		3%	2 %	1 %
yeasticidal (C.albicans)	EN 13624	dirty conditions		1 %	0,5 %	
	EN 14562	dirty conditions		2 %	1 %	
active against enveloped viruses	EN 14476	clean conditions		1 %		
	EN 14476	dirty conditions		3 %		
	EN 17111	dirty conditions		1 %		
tuberculocidal (M. terrae)	EN 14348	clean conditions				5 %
APPLICATION RECOMMENDATION FOR DISINFECTING CLEANING IN THE ULTRASONIC BATH						
bactericidal, yeasticidal	DGHM/VAH	dirty conditions	1 %			

1 – including phase 2 stage 1 and phase 2 stage 2 tests (quantitative suspension tests and practical germ carrier tests)

2 – EN 13624, EN 13727, EN 14561, EN 14562, VAH M 8, VAH M 15

INSTRUMENTS / MANUAL REPROCESSING

PERFEKTAN® TB



PREVENT

PURCHASING INFORMATION

4 / 4

Product	Single unit	Delivery unit	Content	REF
PERFEKTAN TB	Dosage bottle	10	1 L	00-122-010E
	Bottle	6	2 L	00-122-020E
	Canister	3	5 L	00-122-050E

National information may differ. For further information, please contact our subsidiary or your local dealer.

The availability of the products and container sizes depend on a completed national registration.

pH 0 1 2 3 4 5 6 7 8 9 10 11 12 13 14  CE 0123

CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:

<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.