

THERMO CLEAR



CLEAR-RINSING AGENT FOR THERMAL BEDPAN CLEANING MACHINES

1 / 3

Prevents calcification of bedpan washing machines
Ensures spotless drying of the bedpans
Low foaming



INSTRUMENTS / AUTOMATED REPROCESSING

THERMO CLEAR



PRODUCT DESCRIPTION

2 / 3

THERMO CLEAR is suitable for acid thermal reprocessing of bedpans in thermal bedpanwashers. It ensures spotless drying of the bedpan, prevents the cleaning machine for bedpans from calcification and is low-foaming. **THERMO CLEAR** is biodegradable acc. to EU detergent regulation because it is based on citric acid.

APPLICATIONS AND NOTES

Range of application

THERMO CLEAR is suitable for use in thermal bedpan washers of common manufacturers. Do not use with non-thermal and non acid-resistant bedpan washers.

Dosage

0.1 - 0.3 % (1 - 3 ml/L), depending on machine type and water hardness (see instructions of manufacturer). Do not start addition below 30 °C. Dosage is to be done by qualified personnel.

Application

Place materials into the machine according to instructions (e.g. hygiene plan). Every surface and cavity must be accessible for the washing liquor. Use THERMO CLEAR for post-rinsing according to the instructions of the bedpan manufacturer. Attention note for program setting: Do not add at temperatures below 30 °C.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Do not use with non-thermal and non-acid-resistant cleaning machines for bedpans. THERMO CLEAR must not be mixed with other process chemicals.

Composition

Organic acids < 30 %, complexing agents, corrosion inhibitors.

Material compatibility

The product is compatible with all common metals and plastics used in bedpan preparation, such as stainless steel or polypropylene (PP).

Precautionary and hazard statements

Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Dispose of contents/container to approved disposal company or local collection.

For professional use only by personnel with corresponding specialist knowledge according to national directives.

Read label before use.

INSTRUMENTS / AUTOMATED REPROCESSING

THERMO CLEAR



PURCHASING INFORMATION

3 / 3

Product	Single unit	Delivery Unit	Content	REF
THERMO CLEAR	Canister	3	5 L	00-127-050AT
	Flat canister	3	5 L	00-127-050

National information may differ. For further information, please contact our subsidiary or your local dealer.
The availability of the products and container sizes depend on a completed national registration.



CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:
<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.