

DESCOPRENT



PREVENT

INSTRUMENT DISINFECTION FOR DENTAL IMPRESSION MATERIALS

1 / 4

No influence on dimensional accuracy
Ready to use and gentle on materials
Good wetting properties, easily rinsed off



DENTAL

DESCOPRENT



PREVENT

PRODUCT DESCRIPTION

2 / 4

Ready-to-use solution for fast and lasting immersion disinfection of impressions and dental workpieces based on QAV and aldehydes. **DESCOPRENT** disinfects

without negatively affecting the dimensional behavior and surface quality. It behaves color- and odor-neutral and is formaldehyde- and phenol-free.

APPLICATIONS AND NOTES

Range of application

Instrument disinfection for dental impression materials.

Application

For professional use only. Thoroughly rinse the impression material or dental technical workpiece under running water of drinking water quality immediately after removing it from the mouth. Immerse items made from alginates and hydro-colloids shortly in disinfection solution and let dry for 15 minutes. Immerse all other materials into the disinfection solution (observe contact times). Rinse with demineralised water afterwards. The manufacturer's instructions for the impression material must be observed.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Change the working solution immediately if it is visibly soiled, but at least once every working day. DESCOPRENT must not be mixed with other process chemicals.

Composition

100 g contain: 0.48 g dimethylbenzylalkylammonium chloride, 0.2 g ethanedial, 0.01 g pentanedial

Precautionary and hazard statements

Contains Glyoxal, Glutaraldehyde. May produce an allergic reaction. Safety data sheet available on request.

Read label before use.



SPECTRUM OF EFFICACY AND CONTACT TIMES

3 / 4

SPECTRUM OF EFFICACY AND CONTACT TIMES			1 min	5 min	15 min
application recommendations for instrument disinfection					
bactericidal, yeasticidal ¹	VAH/EN ²	dirty conditions		•	
active against enveloped viruses	EN 14476 EN 17111	dirty conditions			•
additional test results					
bactericidal	EN 13727	dirty conditions	•		
	EN 14561	dirty conditions		•	
yeasticidal (C. albicans)	EN 13624	dirty conditions	•		
	EN 14562	dirty conditions		•	
active against enveloped viruses	EN 17111	clean and dirty condition			•

¹ – including phase 1 step 1 and phase 2 step 2 tests (quantitative suspensions tests and practical carrier tests)

² – EN 13624, EN 13727, EN 14561, EN 14562

DESCOPRENT



PURCHASING INFORMATION

4 / 4

Product	Single unit	Delivery unit	Content	REF
DESCOPRENT	Canister	3	5 L	00-813-050 00-813-050EXD

National information may differ. For further information, please contact our subsidiary or your local dealer.
The availability of the products and container sizes depend on a completed national registration.



CE 0123

CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:
<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001,
has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.